

BIOLUX_{PRO} DENT



Post Braces Adjustment

1. Connect your **DENT** device to your phone or to a power bank.
2. Place in your mouth with the tab facing up.
3. Use it for 5 minutes, three times a day.
4. Repeat for 48 hours after adjusting your braces.
5. When the indicated time finishes, disconnect and remove from your mouth.
6. Wash your device with plain water and soap or toothpaste, without wetting the connection indicated with an arrow in this protocol.



Proper hydration enhances results
This does not suggest or replace any medical treatment

