





BIOLUX_{PRO}

MAX WITH TENS



Thigh

1. Place BioluxPro over the thigh, adjusting it with the elastic if necessary.
2. You can place the electrodes on the muscle(s) you want to relax or stimulate. This is complementary and not mandatory to achieve the expected results from the red light.
3. Press  select ALL.
4. Press  select 20 mins.
5. Press  select L10.
6. If you placed the electrodes, press **M** for 3 seconds to turn them on, choose the program that suits you best and use +/- to increase or decrease the intensity.
7. If you feel too much heat during the process, press **F** at 40Hz. to pulse the light.
8. Repeat this process two or three times a day for three days.



Proper hydration enhances results

This does not suggest or replace any medical treatment

