

# BIOLUX<sup>PRO</sup>

## LENS



### Tired Eyesight

1. Wash your face before using, it must be cream, sun block, makeup and dust free. Leave your cheekbone area a little damp.
2. Press **⏻/Mode** three seconds to turn your **LENS** on.
3. Press **+/-** to turn the three lights on.
4. Place the **LENS** below your eyes, touching your skin. Remember your cheekbones are damp.
5. Leave it on for 5 minutes, use them 3 times a day.
6. Turn them off pressing **⏻/Mode**
7. There might be a slight redness in the area for a few minutes.



Proper hydration enhances results

This does not suggest or replace any medical treatment

