



BIOLUX PRO MINI



Wisdom Tooth

1. Place the band around the mouth, adjusting it with the elastic at the back of the neck.
2. Press  twice to activate the red and infrared light.
3. Allow the device to turn off (15 minutes), and repeat this procedure three times a day (morning, afternoon, and night) for 72 hours. Additionally, at night, add 5 minutes by pressing  once with the mouth open.
4. If discomfort persists, consult a healthcare professional.



Proper hydration enhances results

This does not suggest or replace any medical treatment

